ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human CD138, Clone MI15, is intended for use in immunohistochemistry (IHC). This antibody is a useful aid for the classification of multiple myeloma (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Synonym for antigen	Syndecan-1 (1).
Summary and explanation	<p>CD138 is a transmembrane proteoglycan, bearing heparan sulfates and sometimes chondroitin sulfates. It consists of a 30.5 kDa core protein and five glycosaminoglycans. CD138 is a member of the syndecan family. The reported Mr of CD138 ranges from 85 000 to 92 000, presumably depending on the glycosaminoglycan pattern of the molecule (1).</p> <p>Within the hematopoietic system, CD138 is mainly confined to late stages of B-cell differentiation. It is typically expressed at high cellular density in normal and malignant plasma cells (4, 5) and in lymphoplasmacytoid cells (4). In mature tissue CD138 is expressed by simple and stratified epithelia, fibroblasts, stratified keratinocytes and endothelial cells (1).</p> <p>The function of CD138 is largely unknown. However, the intracellular tail of CD138 seems to bind together the cytoskeleton with cellular components (1), and the extracellular part of the molecule seems to bind different ligands. Presumably, CD138 plays a role in cellular functions, such as proliferation, programmed cell death, cell-matrix and cell-cell adhesion (5), for example it mediates adhesion of myeloma cells to type I collagen (6).</p> <p>CD138 expression may be reduced during malignant transformation of various epithelia (7), and may be lost (shed) by myeloma cells entering into apoptosis (1, 6). CD138 may be lost upon malignization of several stratified epithelia (7).</p> <p>Antibodies to CD138 may also be useful for the subclassification of diffuse large B-cell lymphoma (3).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN₃.</p> <p><u>Clone:</u> MI15. <u>Isotype:</u> IgG1, kappa.</p> <p><u>Mouse IgG concentration:</u> See label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	A mixture of U266 and XG-1 human myeloma cell lines (2, 8).
Specificity	<p>Anti-CD138, MI15, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Kobe 1996), and studies by a number of different laboratories confirmed its reactivity with CD138 (1).</p> <p>The epitope recognized by anti-CD138, MI15, is found within the ectodomain (extracellular part) of the CD138 core protein. Blocking experiments have revealed that the epitope recognized by MI15 seems to be partly overlapping the epitope recognized by another antibody to CD138, clone B-B4 (1, 5). Anti-CD138, MI15, reacts strongly with CD138 expressed by multiple myeloma-derived cell lines, such as U266 and XG-1 (8).</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used on formalin-fixed, paraffin-embedded tissue sections. Antigen retrieval of deparaffinized tissue sections by heating in Dako Target Retrieval Solution, Code S1700 or S1699 is mandatory.</p> <p>The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human CD138, Code M7228, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections using 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931 diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen.</p> <p><u>Visualization:</u> Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p>

Limitations	Owing to the sensitivity of CD138 towards proteolytic cleavage, the processing of samples, such as freezing or storage before staining with the antibody, may play an important role in the outcome of the staining (1).
Performance characteristics	Cells labelled by the antibody display staining confined to the cell surface membrane. <u>In blood</u> : No significant reaction is seen with leucocytes from normal peripheral human blood (< 5%). <u>In bone marrow</u> : No significant reaction is seen in normal bone marrow (< 5%). CD138 is absent from CD34 positive stem cells (1). In bone marrow from multiple myeloma patients, all plasma cell types are labelled, including reticular, polymorphous, asynchronous and blastic plasma cells (2).

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

1. Wijdenes J, Clément C, Klein B, Dore J-M. New B-cell CD antigens. BC29: CD138 (syndecan-1) workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 249-52.
2. Costes V, Magen V, Legouffe E, Durand L, Baldet P, Rossi J-F, et al. The M115 monoclonal antibody (anti-syndecan-1) is a reliable marker for quantifying plasma cells in paraffin-embedded bone marrow biopsy specimens. Hum Pathol 1999;30:1405-11.
3. Colomo L, López-Guillermo A, Perales M, Rives S, Martínez A, Bosch F, et al. Clinical impact of the differentiation profile assessed by immunophenotyping in patients with diffuse large B-cell lymphoma. Blood 2003;101:78-84.
4. Sebestyén A, Berczi L, Mihalik R, Paku S, Matolcsy A, Kopper L. Syndecan-1 (CD138) expression in human non-Hodgkin lymphomas. Br J Haematol 1999;104:412-9.
5. Gattei V, Godeas C, Degan M, Rossi FM, Aldinucci D, Pinto A. Characterization of anti-CD138 monoclonal antibodies as tools for investigating the molecular polymorphism of syndecan-1 in human lymphoma cells. Br J Haematol 1999;104:152-62.
6. Inki P, Jalkanen M. The role of syndecan-1 in malignancies [review]. Ann Med 1996;28:63-7.
7. Horvathova M, Gaillard JP, Liautard J, Duperray C, Lavabre-Bertrand T, Bourquard P, et al. B30.6. Identification of novel and specific antigens of human plasma cells by mAb. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 713-4.

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Agilent Technologies Singapore (International) Pte Ltd.