**Monoclonal Mouse
Anti-Human CD68**
Clone KP1
Code M0814

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human CD68, Clone KP1, is intended for use in immunohistochemistry (IHC). The antibody labels macrophages and other members of the mononuclear phagocyte lineage and is a useful aid for the classification of neoplasms of myeloid and macrophage/monocyte origin (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and Explanation	<p>CD68 is a highly glycosylated lysosomal membrane protein with an Mr of 110 000. The CD68 protein belongs to a family of lysosomal glycoprotein (LGP)/plasma membrane shuttling proteins that play a role in endocytosis and/or lysosomal trafficking. CD68 is expressed strongly in cytoplasmic granules, and weakly on the surface of macrophages, monocytes, neutrophils, basophils and NK-cells. Additionally, CD68 is expressed by some peripheral blood B cells and may be weakly expressed in B-cell type acute lymphoblastic leukemia (B-ALL) cells. CD68 can also be found in the cytoplasm of non-hematopoietic tissues, especially the liver, and renal glomeruli, and tubules (2). Unlike many other CD leucocyte antigens, the CD68 molecule is antigenically very heterogenous, and different antibodies to CD68 show different cellular reactivities (3).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L Na₃N.</p> <p><u>Clone</u>: KP1 (4). <u>Isotype</u>: IgG1, kappa.</p> <p><u>Mouse IgG concentration</u>: see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Lysosomal fraction of human lung macrophages (4).
Specificity	<p>The antibody was clustered as anti-CD68 at the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Vienna in 1989 (5).</p> <p>SDS-PAGE analysis of immunoprecipitates formed between the antibody and ¹²⁵I-labeled lysates from human spleen with B-cell lymphoma rich in macrophages shows reaction with a 110 kDa polypeptide, corresponding to CD68 (4).</p> <p>In Western blotting of extracts of lung, spleen and U937 cells, diffuse 110, 70 and 40 kDa bands were detected when using reducing conditions. Under non-reducing conditions the spleen extract showed an additional 220 kDa band (4).</p>
Precautions	<ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (Na₃N), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<p><u>Paraffin sections</u>: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formal saline, B5 (4) or formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found less efficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution</u>: Monoclonal Mouse Anti-Human CD68, Code M0814, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control</u>: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p> <p><u>Visualization</u>: Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p>
Product specific limitations	The antibody shows a lack of specificity in frozen sections where it labels a number of normal and neoplastic structures of lymphoid and epithelial components (6). Melanomas are commonly immunoreactive for histiocytic markers and the KP1 antibody labels about 68% of melanomas, typically providing diffuse and weak staining. Accordingly, the antibody should be used cautiously (7).

Staining interpretation	Myelomonocytic cells labeled by the antibody display a diffuse or granular cytoplasmic staining.
Performance characteristics	<p><u>Normal tissues:</u> In normal peripheral blood smears all monocytes and most granulocytes are labeled by the antibody. Tissue macrophages in a wide range of paraffin-embedded tissues are labeled, including lung, germinal center and bone marrow macrophages, and Kupffer cells in the liver. In the bone marrow also myeloid precursors, and many mature granulocytes are strongly labeled (4). Further, microglial cells in the brain (8), osteoclasts in bone, kidney glomeruli and mast cells are strongly labeled, while moderate reaction is observed with kidney tubules, hepatocytes and occasional lymphoid cells (6). Langerhans' cells, interdigitating reticulum cells, and follicular dendritic cells (FDC) are not labeled, except for weak labeling observed in a minority of FDCs in dermatopathic lymphadenopathy (4).</p> <p><u>Abnormal tissues:</u> Acute myeloid leukemias of M1 to M5 type are strongly labeled by the antibody (6). In another study, 20/20 neoplasms of myeloid, myelomonocytic and presumed macrophage derivation showed strong and extensive cytoplasmic reactivity with the antibody. 14/41 B-lineage lymphomas and leukemias were also labeled, but the labeling was usually confined to small dots. The labeled B-cell neoplasms were almost all small-cell proliferations (1). In 23/36 (64%) plasma cell hyperplasias, more than 1% of the plasma cells were labeled by the antibody (9). Of 43 primary and metastatic melanomas, 86% were labeled weakly by the antibody (7). All of 22 T-cell lymphomas were not labeled, as also 12/12 CD30+ anaplastic large-cell lymphomas (1).</p>

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
References/ Bibliographie/ Literatur

1. Warnke RA, Pulford KAF, Pallesen G, Ralfkiaer E, Brown DC, Gatter KC, et al. Diagnosis of myelomonocytic and macrophage neoplasms in routinely processed tissue biopsies with monoclonal antibody KP1. Am J Pathol 1989;135:1089-95.
2. Goyert SM. MC12. CD68 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1015-16.
3. Falini B, Flenghi L, Pileri S, Gambacorta M, Bigerna B, Durkop H, et al. PG-M1: A new monoclonal antibody directed against a fixative-resistant epitope on the macrophage-restricted form of the CD68 molecule. Am J Pathol 1993;142:1359-72.
4. Pulford KAF, Rigney EM, Micklem KJ, Jones M, Stross WP, Gatter KC, et al. KP1: a new monoclonal antibody that detects a monocyte/macrophage associated antigen in routinely processed tissue sections. J Clin Pathol 1989;42:414-21.
5. Micklem K, Cordell J, Rigney E, Simmons D, Pulford K, Stross P, et al. M13.1. A macrophage-associated antigen defined by five mAb. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 843-46.
6. Cordell JL, Falini B, Flenghi L, Jones DB, Pileri S, Radzun HJ, et al. M15. CD68 cluster workshop report. In: Schlossman SF, Bousmell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p.925-7.
7. Pernick NL, DaSilva M, Gangi MD, Crissman J, Adsay V. "Histiocytic markers" in melanoma. Mod Pathol 1999;12:1072-7.
8. Aoki T, Kobayashi K, Isaki K. Microglial and astrocytic change in brains of Creutzfeldt-Jakob disease: an immunocytochemical and quantitative study. Clin Neuropathol 1999;18:51-60.
9. Beschorner R, Horny H-P, Petrucci UR, Kaiserling E. Frequent expression of haemopoietic and non-haemopoietic antigens by reactive plasma cells: an immunohistochemical study using formalin-fixed, paraffin-embedded tissue. Histol Histopathol 1999;14:805-12.

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Agilent Technologies Singapore (International) Pte Ltd.