**Monoclonal Mouse
Anti-Human
CD79 α
Clone JCB117
Code M7050**

ENGLISH

Intended use

For in vitro diagnostic use.

Monoclonal Mouse Anti-Human CD79 α , Clone JCB117, is intended for use in immunohistochemistry (IHC). The antibody labels B cells and is a useful aid for the classification of B-cell neoplasms of all maturation stages (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

Synonyms for antigen

Ig- α , mb-1 (2).

Summary and explanation

CD79 is a disulphide-linked transmembrane heterodimer with a molecular mass of 82-95 kDa belonging to the immunoglobulin superfamily. It comprises the two glycoproteins, CD79 α , with a molecular mass of 40-45 kDa, and CD79 β , with a molecular mass of 37 kDa. CD79 is non-covalently associated with surface Ig, forming the B-cell receptor complex, which is required for antigen recognition. CD79 is essential for efficient surface expression of the B-cell receptor, and is also necessary for signal transmission into the cytoplasm following antigen-binding to surface Ig. The cytoplasmic portion of CD79 α and CD79 β contains several kinases. Phosphorylation signals via these kinases result in B-cell activation, differentiation and, in some cases, apoptosis (2-4).

The expression of CD79 is largely restricted to B-lineage cells, but CD79 α may be coexpressed with CD3 in a proportion of T-lymphoblastic leukemia/lymphoma (5). In precursor B cells, the CD79 protein chains are already expressed in the cytoplasm (CyCD79). Surface expression of CD79 begins at the pro-B cell stage and persists throughout the B-cell differentiation. CD79 α is more strongly expressed by B cells in the follicular mantle zone than by germinal centre B cells, suggesting that activation of mature B cells downregulates CD79 expression. CD79 expression ceases around the onset of plasma cell differentiation, with only a proportion of plasma cells containing CD79 (2, 4).

Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.

Reagent provided

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaNa₃.

Clone: JCB117 (1). IsoType: IgG1, kappa.

Mouse IgG concentration: see label on vial.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Immunogen

Recombinant protein containing part of the extracellular portion of the human CD79 α glycoprotein (1).

Specificity

The antibody was clustered as anti-CD79 α at the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Kobe in 1996 (4).

In Western blotting of Ramos B-cell lysate, the antibody labels a band corresponding to CD79 α under reducing conditions, and a band corresponding to CD79 $\alpha\beta$ under non-reducing conditions (1).

SDS-PAGE analysis of immunoprecipitates formed between the antibody and ¹²⁵I-labeled, mildly denatured, and reduced lysate of Ramos B cells, shows weak reaction with a 43 kDa polypeptide corresponding to CD79 α .

The antibody was tested by flow cytometry against normal peripheral B-cell lines, and against mature and immature B-cell lines (Ramos, Raji, and NALM-1), and was consistently negative, indicating that the epitope recognized is not present, or is not accessible, on living cells. However, the antibody specifically labels B cells in paraffin-embedded tissue sections (4).

Precautions

1. For in vitro diagnostic use.
2. For professional users.
3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin sections: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin, B5 fixative, Bouin's fixative, or acetic formalin (1). Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.

Staining procedure

These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.

Dilution: Monoclonal Mouse Anti-Human CD79 α , Code M7050, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.

Visualization: Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Product-specific limitations

Of 149 cases of T-lymphoblastic leukemia/lymphoma, 10% showed sufficiently extensive CD79 α labeling with the antibody (more than 90% labeled cells) to cause assessment problems. All of the 149 T-cell neoplasms were, however, labeled for the T-cell marker CD3, while none of 68 cases of B-lymphoblastic leukemia/lymphoma were CD3 labeled (5).

Staining interpretation

Cells labeled by the antibody display staining of the cell membrane and/or the cytoplasm.

Performance characteristics

Normal tissues: Normal plasma cells in tissue samples are strongly labeled with the antibody (1). In tonsil the antibody strongly labels B cells in the germinal center and the mantle zone, and B cells present in T-cell areas.

Abnormal tissues: The antibody labeled all of 331 B-cell neoplasms, including 41 lymphoblastic lymphomas/leukemias, 28 small lymphocytic, 36 lymphoplasmacytoid, 17 mantle cell, 53 follicular, 29 MALT, 95 large cell, and 7 Burkitt's lymphomas. Additionally, 15/15 hairy cell leukemias, 13/15 B-cell anaplastic large cell lymphomas, and 10/20 myelomas/plasmacytomas were labeled with the antibody. In the same study all of 98 T-cell and non-lymphoid neoplasms were not labeled with the antibody, including 9 T-cell lymphoblastic lymphomas/leukemias, 10 mycosis fungoides, 32 peripheral T-cell lymphomas, 8 angioblastic T-cell lymphomas, 11 T-cell anaplastic large cell lymphomas, and 28 acute myeloid leukemias (1). Two later studies (5, 6), have shown immunoreaction with T-cell neoplasms. Thus, in 149 cases of CD3+ T (precursor)-acute lymphoblastic leukemia/lymphoma, 14 cases expressed CD79 α on more than 90% of the blast cells as determined with the antibody. In 55 cases 10-90% of the blast cells were labeled (5). Further, in 4/94 of enteropathy-type intestinal T-cell lymphomas and 1/11 nasal NK/T-cell lymphomas which were CD3+, the majority of tumor cells were labeled by the antibody (6).

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Agilent Technologies Singapore (International) Pte Ltd.