**Polyclonal Rabbit
Anti-Human
Calcitonin**
Code A0576

ENGLISH	
Intended use	<p>For in vitro diagnostic use.</p> <p>Polyclonal Rabbit Anti-Human Calcitonin, Code A0576, is intended for use in immunohistochemistry (IHC). The antibody labels calcitonin-producing C-cells. Differential classification of tumors is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>Calcitonin (CT) is a 32 amino acid peptide hormone derived from a precursor consisting of 136 amino acids and an N-terminal leader sequence. CT contains a carboxyterminal proline amide and a disulphide bridge between cysteine residues at position 1 and 7. The sequence is highly conserved in a variety of species within the N-terminal loop region but demonstrates divergence in the rest of the sequence. CT is produced in the parafollicular C-cells of the thyroid and acts through its receptors (CTR) inhibiting osteoclast-mediated bone resorption and increasing calcium excretion by the kidney. Accumulating evidence supports the existence of non-thyroidal CT-like peptides in tissues such as brain, prostate and uterus, which suggest additional physiological roles of CT (1). Medullary thyroid carcinoma (MTC) may be sporadic or familial, the latter has an autosomal dominant pattern of inheritance (2, 3).</p> <p>Antibodies to calcitonin have been shown to aid in the classification of MTC (2).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage, Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mol/L NaN₃.</p> <p><u>Protein concentration:</u> See label on vial.</p> <p>The titre variation between different lots of A0576 is less than 10%. This is achieved by adjusting the titre of each individual lot to match the titre of an antibody reference preparation kept at -80 °C.</p>
Immunogen	Synthetic human calcitonin, aa 1-32, coupled to bovine serum albumin (BSA).
Specificity	<p>The antibody reacts with human calcitonin. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma and bovine serum proteins.</p> <p>The antibody labels human thyroid C-cells and other calcitonin-containing cells strongly. No non-specific staining has been observed.</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Polyclonal Rabbit Anti-Human Calcitonin, Code A0576, may be used at a dilution range of 1:400-1:800 when applied on formalin-fixed, paraffin-embedded sections of human thyroid tissue and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), Code X0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p>

Visualization: Dako EnVision+/*HRP* kits, e.g. Code K4009, are recommended. Follow the procedure enclosed with the selected visualization kit.

Staining interpretation Cells labelled by the antibody display a cytoplasmic staining pattern.

Performance characteristics Normal tissues: The antibody labels parafollicular C-cells of the thyroid gland. No labelling is observed in the adrenal gland (4).

Abnormal tissues: The antibody labels medullary thyroid carcinomas (2). In addition, the antibody labeled 6/28 pheochromocytomas (4), 7/37 carcinoids of the duodenum (5), 2/24 gastric carcinoids (6), and tumor cells in 11/15 cases – and surface cells in 2/15 cases of sclerosing hemangiomas of the lung (7). No labeling was observed in 27 cases of adrenocortical carcinomas (4).

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
References / Bibliographie / Literaturnachweise

1. Pondel M. Calcitonin and calcitonin receptors: bone and beyond. *Int J Exp Pathol* 2000;81:405-22.
2. Hayashida CY, Alves VAF, Kanamura CT, Ezabella CL, Abelin NM, Nicolau W, et al. Immunohistochemistry of medullary thyroid carcinoma and C-cell hyperplasia by an affinity-purified anti-human calcitonin antiserum. *Cancer* 1993;72:1356-63.
3. Giuffrida D, Gharib H. Current diagnosis and management of medullary thyroid carcinoma. *Ann Oncol* 1998;9:695-701.
4. Komminoth P, Roth J, Schröder S, Saremaslani P, Heitz PV. Overlapping expression of immunohistochemical markers and synaptophysin mRNA in pheochromocytomas and adrenocortical carcinomas. Implications for the differential diagnosis of adrenal gland tumors. *Lab Invest* 1995;4:424-31.
5. Burke AP, Federspiel BH, Shekitka KM, Helwig EB. Carcinoids of the duodenum. A histologic and immunohistochemical study of 65 tumors. *Am J Surg Pathol* 1989;13:828-37.
6. Thomas RM, Baybick JH, Elsayed AM, Sobin LH. Gastric carcinoids. An immunohistochemical and clinicopathologic study of 104 patients. *Cancer* 1994;73:2053-8
7. Xu H-M, Li W-H, Hou N, Zhang S-G, Li H-F, Wang S-Q, et al. Neuroendocrine differentiation in 32 cases of so-called sclerosing hemangioma of the lung: Identified by immunohistochemical and ultrastructural study. *Am J Surg Pathol* 1997;9:1013-22.

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Agilent Technologies Singapore (International) Pte Ltd.