**Monoclonal Mouse  
Anti-Human Desmin**  
Clone D33

**Code M0760**

**ENGLISH**

<b>Intended use</b>	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Desmin, Clone D33, is intended for use in immunohistochemistry (IHC). The antibody labels smooth and striated muscle cells as well as mesothelial cells, and is a useful aid for the classification of rhabdomyosarcoma (1, 2), leiomyoma (3, 4), and mesothelioma (5). Differential classification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.
<b>Summary and explanation</b>	Desmin belongs to the class III of intermediate filaments, constituting part of the cytoskeleton, and is the characteristic intermediate filament of all three types of muscle cells (skeletal, cardiac and smooth muscle). Desmin is a 53-kDa protein encoded by nine exons of a gene located at 2q35. Desmin forms a cytoskeletal network across the muscle fibre bordering at the plasma and nuclear membrane and is particularly localized to the subplasmalemmal region and the Z-band (6). <i>Refer to <a href="#">Dako General Instructions for Immunohistochemical Staining</a> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</i>
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> . <u>Clone:</u> D33. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
<b>Immunogen</b>	Desmin purified from human muscle.
<b>Specificity</b>	In Western blotting of cytoskeletal preparations of MCF-7 cells, SCC-4 cells RT-112 cells, skin epidermis, primary culture of renal carcinoma cells, leiomyoma, fetal heart, and crude extract from human white matter of brain, the antibody labels only a ~53 kDa band in the leiomyoma and fetal heart preparations, proving the reactivity with desmin, and demonstrating the absence of reactivity with other types of intermediate filaments tested (7).
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For in vitro diagnostic use.</li> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> <li>Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin or B5-fixative (1). Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections:</u> The antibody can be used for labeling frozen sections (7). The user must validate the staining procedure.
<b>Staining procedure</b>	These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms. <u>Dilution:</u> Monoclonal Mouse Anti-Human Desmin, Code M0760, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human colon and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0 and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809. <u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimen. <u>Visualization:</u> Dako EnVision+HRP kits, e.g. Codes K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.
<b>Staining interpretation</b>	Cells labeled by the antibody display a cytoplasmic staining pattern. The staining may show a fibrillar aspect (7).
<b>Performance characteristics</b>	<u>Normal tissues:</u> The antibody labels vascular smooth muscle cells and esophageal smooth muscle cells (3). In 8/10 fixed tissues, the antibody labeled mesothelial cells (5). Except for cells of muscle origin, normal bone marrow, kidney, breast, prostate gland, brain (gray matter), jejunum, thyroid gland, lung, colon, stomach, ureter and tonsil are negative (1). In fetal frozen tissues, the antibody labels smooth muscle cells in blood vessels and the small intestine, skeletal muscle cells in the tongue, cardiac muscle cells, stromal cells of the medulla

of the kidney, the decidua, placental villi and the umbilical cord, submesothelial stromal cells of the pericardium, as well as mesothelial cells (7).

**Abnormal tissues:** Of human rhabdomyosarcomas, the antibody labeled 4/4 tumors (1). In the human esophagus the antibody labeled 35/35 leiomyomas, 1/3 leiomyosarcomas and 3/17 stromal tumors (3). In addition, the antibody labeled 33/33 highly cellular leiomyomas of the uterus (4). Also, 8/16 malignant, diffuse mesotheliomas of the pleura, were labeled by the antibody (5). In human vaginal and endocervical polyps, the antibody labeled both background stromal cells and the atypical stromal cells within the subepithelial zone (8). The antibody showed no labeling in 10 cases each of neuroblastoma, lymphoma, Wilms' tumor, primitive neuroectodermal tumor, and retinoblastoma (2), nor did it label 4 cases of small-cell carcinoma of lung, 3 cases of Ewing's sarcoma (1), and 6 cases of endometrial stromal nodules (4).

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Agilent Technologies Singapore (International) Pte Ltd.