**Polyclonal Rabbit
Anti-Human
Gastrin
Code A0568**

ENGLISH	
Intended use	<p>For in vitro diagnostic use.</p> <p>Polyclonal Rabbit Anti-Human Gastrin, Code A0568, is intended for use in immunohistochemistry (IHC). This antibody labels gastrin-expressing cells and is a useful aid for the classification of gastrin-secreting neuroendocrine tumors (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>Gastrin is a peptide hormone that is important in the regulation of gastric acid secretion and mucosal cell proliferation (3). Progastrin is synthesized in G cells in the antropyloric mucosa, extraantral G cells and in endocrine cells of the pancreas and the pituitary gland. In G cells the progastrin is concentrated into secretory granules, cleaved in the N- and C-terminal and amidated to form gastrin (4). The two major forms of gastrin are gastrin-34 and gastrin-17 (3, 4). Secreted gastrin acts by binding to the G-coupled CKK_2 receptor and stimulates the release of histamine from enterochromaffin-like cells. The released histamine acts as a paracrine stimulant on parietal cells that release gastric acid into the gut lumen (5).</p> <p>Gastrinomas are gastrin secreting neuroendocrine tumors (6) associated with the Zollinger-Ellison (ZES) syndrome or with the autosomal dominant inherited disorder multiple endocrine neoplasia type 1 (MEN1) syndrome. Sporadic gastrinomas are generally single and most commonly found in the pancreas but can be located to the duodenum. MEN1 associated gastrinomas are most often multiple and located to the duodenum (1, 7).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaN_3.</p> <p><u>Protein concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	Synthetic human gastrin-17 non-sulphated form, coupled to bovine serum albumin.
Specificity	<p>The antibody reacts with human non-sulphated and sulphated forms of gastrin-17 as well as gastrin-34, and cross-reacts about 20% with cholecystokinin octapeptide (CCK-8). Traces of contaminating antibodies have been removed by solid-phase absorption with bovine plasma proteins.</p> <p>The specificity of the antibody has been ascertained as follows:</p> <p><u>Crossed immunoelectrophoresis:</u> No precipitation arch appears when using 12.5 μL A0568 per cm^2 gel area against 2 μL of human plasma or 2 μL bovine serum. Staining: Coomassie Brilliant Blue.</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pretreatment of deparaffinized tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Polyclonal Rabbit Anti-Human Gastrin, Code A0568, may be used at a dilution range of 1:400-1:800 when applied on formalin-fixed, paraffin-embedded sections of stomach and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), Code X0936, diluted to SSA0568CEEFG_02 p. 1/5</p>

the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in Dako Antibody Diluent, Code S0809.

Visualization: Dako EnVision+HRP kits, e.g. Code K4009, are recommended. Follow the procedure enclosed with the selected visualization kit.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.

Cells labeled by the antibody display cytoplasmic staining.

Staining interpretation

Performance characteristics

Normal tissues: In stomach, G-cells in the pyloric antrum show moderate to strong staining reaction. The antibody does not label epithelial cells.

Abnormal tissues: The antibody labeled 12/12 duodenal gastrinomas, 5/5 pancreatic gastrinomas, 8/8 gastrinomas in peripancreatic and periduodenal lymph nodes and liver without known primary tumor (1). In another study, the antibody labeled 3/11 pancreatic endocrine tumors (2).

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
1. Perkins PL, McLeod MK, Jin L, Fukuuchi A, Cho KJ, Thompson NW, et al. Analysis of gastrinomas by immunohistochemistry and in situ hybridization histochemistry. *Diag Mol Path* 1992;1:155-4.
2. Wang G, Yao J, Worah S, White J, Luna R, Wu TT et al. Comparison of genetic alterations in neuroendocrine tumors: frequent loss of chromosome 18 in ileal carcinoid tumors. *Modern Pathol* 2005;18:1079-87.
3. Thomas RP, Hellmich MR, Townsend CM and Evers BM. Role of gastrointestinal hormones in the proliferation of normal and neoplastic tissues. *Endocrine Rev* 2003;24:571-99.
4. Dickinson CJ. Relationship of gastrin processing to colon cancer. *Gastroenterology* 1995;109:1384-8
5. Dimaline R, Varro A. Attack and defence in the gastric epithelium-a delicate balance. *Exp Physiol* 2007;92:591-601.
6. de Herder W. Biochemistry of neuroendocrine tumors. *Best Pract Res Clin Endocrinol Metab* 2007;21:33-41.
7. Pritchard DM. Pathogenesis of gastrinomas associated with multiple endocrine neoplasia type 1. *Gut* 2007;56:606-7.

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Agilent Technologies Singapore (International) Pte Ltd.