**Polyclonal Rabbit  
Anti-Human  
Kappa Light Chains  
Code No. / Code/ Code-Nr. A 0191**

## ENGLISH

### Intended use

For in vitro diagnostic use.

Polyclonal Rabbit Anti-Human Kappa Light Chains, code No. A 0191, is intended for use in immunocyto-chemistry. The antibody is also well-suited for gel immunoprecipitation techniques, including immunofixation and immunoblotting (1).

In immunocytochemistry, the antibody labels plasma cells and related lymphoid cells containing kappa light chains, and it is a useful tool for the classification of patients with monoclonal gammopathies (2) and amyloidosis (3, 4). Additionally, the antibody may be used for distinguishing neoplastic monoclonal proliferation from reactive hyperplasia of B cells (5-7). Differential identification is aided by the results from a panel of antibodies.

Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

### Introduction

The human immunoglobulins basically consist of two identical heavy chains (Mr about 50 000) and two identical light chains (Mr about 20 000). The four chains are covalently linked together by disulphide bonds. The light chains are either kappa or lambda. The five immunoglobulins, IgA, IgD, IgE, IgG and IgM, differ regarding the heavy chains, and IgA and IgM also occur in polymeric forms (8).

Individual B cells express either kappa or lambda light chains, but never both. In a polyclonal population the ratio of kappa to lambda bearing B cells is 2:1, and the occurrence of a mixture of kappa and lambda light chain-bearing cell types suggests polyclonality and a reactive or non-neoplastic proliferation of B cells (9). Diseases, such as multiple myeloma and B-cell lymphoma, are characterized by the proliferation of monoclonal neoplastic plasma cells producing only one type of light chain. Thus, the demonstration of a kappa or lambda light chain-restricted cell population is very useful in the histopathological diagnosis of these diseases (4). Primary AL (amyloid light-chain) amyloidosis is also a plasma cell disorder where the demonstration of monotypic kappa or lambda light chain in the amyloid deposits is diagnostically important and separates this disease from AA amyloidosis where immunoglobulin light chain deposits do not occur (3).

### Reagent provided

Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaNa<sub>3</sub>.

**Protein concentration g/L:** See label on vial. **Antibody titre (SRI):** 600 mg/L (10).

As no international kappa reference preparation is available, a polyclonal kappa preparation (isolated from human IgG) has been used for the titre determination. The titre variation between different lots of A 0191 is less than 10%. This is achieved by adjusting the titre of each individual lot to match the titre of an antibody reference preparation kept at -80 °C.

### Immunogen

Polyclonal immunoglobulin light chains of kappa type isolated from a pool of human sera.

### Specificity

The antibody reacts with free kappa chains as well as kappa chains in intact immunoglobulin molecules. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins.

Specificity tests; staining with Coomassie Brilliant Blue:

**Double immunodiffusion:** No reaction with lambda or gamma chains is seen when using 15 µL of antibody against 15 µL of lambda IgG.

**Crossed immunoelectrophoresis:** Only kappa-related precipitates appear when using 12.5 µL antibody per square cm gel area against 2 µL human plasma or 2 µL of a pool of concentrated urines from patients with tubular proteinuria.

Cross-reaction with immunoglobulins from other species may occur. As demonstrated by immunocytochemistry, the antibody cross-reacts with immunoglobulins in dog, horse, pig and rat.

### Precautions

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. The product may be used in different techniques and in combination with different sample types and materials, therefore each individual laboratory should validate the test system applied.

### Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

## IMMUNOCYTOCHEMISTRY

### Specimen preparation

**Paraffin sections:** The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or B5 (2). Pre-treatment of tissues with proteinase K or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, code No. S 1700, Dako Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.

**Frozen sections and cell preparations:** The antibody can be used for labelling acetone-fixed, frozen sections (6).

### Staining procedure

**Dilution:** Polyclonal Rabbit Anti-Human Kappa Light Chains, code No. A 0191, may be used at a dilution range of 1:1000-1:2000 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), code No. X 0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in Dako Antibody Diluent, code No. S 0809.

**Visualization:** DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4008 and K 4010, are recommended. Follow the procedure enclosed with the selected visualization kit.

**Automation:** The antibody is well-suited for immunocytochemical staining using automated platforms, such as the Dako Autostainer

### Performance characteristics

Cells labelled by the antibody display staining of the cell membrane and/or cytoplasm.

**Normal tissues:** In normal tonsil, the antibody labels plasma cells, follicle centre cells and mantle zone cells providing a polyclonal pattern with a clear distinction between positive and negative cells (4).

**Abnormal tissues:** In reactive lymph nodes, a characteristic staining pattern was demonstrated with Dako Anti-Kappa, and an indistinguishable pattern was obtained with Dako Anti-Lambda clearly demonstrating the polyclonality of the B cells in this disorder (6). In plasmacytomas, the abnormal plasma cells were monotypic and either kappa or lambda-positive (5). Using frozen sections, all of 15 nodular lymphomas showed a striking predominance of one light chain type when labelled with Dako Anti-Kappa or Dako Anti-Lambda. Of 31 diffuse lymphomas, 23 had a majority of cells with monotypic light chains, in most cases similar to the nodular lymphomas. All of 16 diffuse large cell lymphomas also showed monotypic light chain labelling (6). In 113 cases of formalin-fixed B-cell non-Hodgkin lymphoma, including several small core biopsy specimens with extremely limited tissue, monotypic light chain expression was demonstrated in 91 (81%) of the cases with a 100% specificity using Dako Anti-Kappa and Anti-Lambda (7). In bone marrow aspirates from 66 patients with serum monoclonal gammopathies, the antibody was used to identify abnormal light chain ratios in 24 cases of gammopathy of undetermined significance, 23 cases of multiple myeloma, 5 cases of amyloidosis, 9 cases of macroglobulinaemia, 1 case of extramedullary plasmacytoma, 3 cases of renal disease, and 1 case of lymphoma. The predominant light chain labelled by the antibody in the biopsy correlated with the monoclonal light chain identified in the serum (2). In 23 cases of light chain deposition disease (LCDD), the antibody labelled renal specimens monotypically in 17 cases, while 5 cases were positive for lambda chains only (3).

## IMMUNOFIXATION

### Dilution guidelines

Antibody: Undiluted, 80 µL is usually sufficient. Sample: Human serum 1+3 or urine (protein concentration about 1 g/L). Sample volume: 5 µL. Plasma samples should not be used as the fibrinogen band may interfere with the interpretation.

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
- Withold W, Reinauer H. An immunoblotting procedure following agarose gel electrophoresis for detection of Bence Jones proteinuria compared with immunofixation and quantitative light chain determination. Eur J Clin Chem Clin Biochem 1995;33:135-8.
- Peterson LC, Brown BA, Crosson JT, Mladenovic J. Application of the immunoperoxidase technic to bone marrow trephine biopsies in the classification of patients with monoclonal gammopathies. Am J Clin Pathol 1986;85:688-93.
- Strøm EH, Fogazzi GB, Banfi G, Pozzi C, Mihatsch MJ. Light chain deposition disease of the kidney. Morphological aspects in 24 patients. Virchows Arch 1994;425:271-80.
- Jessup E. Antigen retrieval techniques for the demonstration of immunoglobulin light chains in formalin-fixed, paraffin-wax embedded sections. UK NEQAS immunocytochemistry news. Winter 1994/95; Issue 4:12-6.
- Taylor CR, Burns J. The demonstration of plasma cells and other immunoglobulin-containing cells in formalin-fixed, paraffin-embedded tissues using peroxidase-labelled antibody. J Clin Pathol 1974;27:14-20.
- Harris NL, Poppema S, Data RE. Demonstration of immunoglobulin in malignant lymphomas. Use of an immunoperoxidase technic on frozen sections. Am J Clin Pathol 1982;78:14-21.
- Marshall-Taylor CE, Cartun RW, Mandich D, DiGuseppe JA. Immunohistochemical detection of immunoglobulin light chain expression in B-cell non-Hodgkin lymphomas using formalin-fixed, paraffin-embedded tissues and a heat-induced epitope retrieval technique. Appl Immunohistochem Mol Morphol 2002;10:258-62.
- Klein J, Hořejši V. Immunology. 2nd ed. Abingdon (UK): Blackwell Science Ltd; 1999. p. 226-7.
- Leong ASY, Cooper K, Leong FJWA. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 218.
- Becker W. Determination of antisera titres using the single radial immunodiffusion method. Immunochem 1969;6:539-46.

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**Spezifitätstests:** Anfärben mit Coomassie © Brillantblau:
Agilent Technologies Singapore (International) Pte Ltd.