**Monoclonal Mouse
Anti-Human
Prostatic Acid Phosphatase
Clone PASE/4LJ**

Code M0792

ENGLISH

| | |
|--------------------------------|---|
| Intended use | <p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Prostatic Acid Phosphatase, Clone PASE/4LJ, is intended for use in immunohistochemistry (IHC). The antibody labels glandular epithelium of the prostate and is a useful aid for the classification of primary and metastatic carcinoma of the prostate. Occasional samples of carcinoma of the bladder, rectal carcinoids, rare cases of other tumors, and macrophages may also be labeled by the antibody (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p> |
| Synonyms for antigen | PAP (1, 2), PSAP (3). |
| Summary and explanation | <p>Prostatic acid phosphatase (PAP) belongs to the group of acid phosphatase isoenzymes (EC 3.1.3.2). It is found in large amounts in the prostate, and the concentration in seminal fluid is approximately 1 g/L (1, 2, 4). PAP is synthesized in the secretory epithelial cells in the normal prostate gland. The physiological function of PAP is not known. Chemically it hydrolyzes a wide variety of phosphomonoesters under acidic conditions and is inhibited by L-tartrate. These features make PAP belong to the group of tartrate-inhibitable acid phosphatases, which also includes lysosomal acid phosphatase (4). Immunohistochemical data suggests that detection of PAP is a relatively useful test for the prostatic origin of an otherwise unclassifiable primary or metastatic neoplasm (3).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p> |
| Reagent provided | <p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L Na₂S₂O₃.</p> <p><u>Clone</u>: PASE/4LJ (1). <u>Isotype</u>: IgG1, kappa.</p> <p><u>Mouse IgG concentration</u>: see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p> |
| Immunogen | Purified prostatic acid phosphatase from human seminal plasma (1). |
| Specificity | <p>In Western blotting of whole seminal plasma, the antibody labels a single band of 52 kDa corresponding to prostatic acid phosphatase (1).</p> <p>In ELISA the antibody specifically reacts with human semen excluding all other bodily excretions, including vaginal fluid, blood, saliva, female urine, nasal discharge, earwax, sweat and faeces.</p> |
| Precautions | <ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (Na₂N₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations. |
| Storage | Store at 2-8 °. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support. |
| Specimen preparation | <p><u>Paraffin sections</u>: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin, formal acetic acid, Bouin's formal sublimate or methacarn. Additionally, bone trephines decalcified in either EDTA or in formic/citric acid may be used (2). Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is recommended. For tissues fixed in formalin, optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations</u>: The antibody can be used for labeling acetone-fixed, frozen sections (2, 6). The user must validate the staining procedure.</p> |
| Staining procedure | <p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution</u>: Monoclonal Mouse Anti-Human Prostatic Acid Phosphatase, Code M0792, may be used at a dilution range of 1:200-1:400 when applied on formalin-fixed, paraffin-embedded sections of human prostate and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control</u>: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p> |

Visualization: Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.

Staining interpretation

Cells labeled by the antibody display a cytoplasmic staining pattern.

Performance characteristics

Normal tissues: In normal prostate, the antibody labels glandular epithelial cells. The antibody also labels the prostatic fibromuscular stroma to a variably degree, though smooth muscle elsewhere is not labeled. Staining is also observed in scattered cells within the urothelium of adult, fetal prostatic and penile urethra. Additionally the antibody labels the loops of Henle, maculae densae and distal tubules in 10/10 fetal and 2/9 adult kidneys. Islet cells in 4/5 specimens of normal pancreas and 4/9 blocks of normal pancreas surrounding a pancreatic tumor were also labeled by the antibody (2).

Apart from the labeling mentioned above, no labeling was observed in the genitourinary system, lympho-proliferative system, musculoskeletal system, cardiovascular and respiratory system, endocrine and glandular tissue and the nervous system. No labeling was seen in 13 fetal tissues and in first, second and third trimester placenta and umbilical cord (2).

Abnormal tissues: The antibody labeled all of 12 benign, 32 malignant, and 11 metastatic prostates. Of 402 non-prostatic neoplasms, 9 were labeled by the antibody, i.e. 6/36 carcinoids, 1/55 ovarian adenocarcinomas, 1/1 carcinosarcoma of the lung and 1/9 islet cell tumors. In the remaining 393 tumors no labeling was observed, thus 22 tumors of the urinary tract, 78 in the genital tract, 62 in endocrine organs, 30 carcinoids, 38 tumors in bone, muscle and soft tissue, 32 in breast, 44 in the gastrointestinal tract, 36 in the respiratory tract, 17 in the lympo-proliferative system, 16 in the skin and 18 in the nervous system were all non-immunoreactive (2).

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
References/ Bibliographie/ Literaturnachweise

- Haines AMR, Larkin SE, Heyderman E. A new monoclonal antibody to human prostatic acid phosphatase suitable for immunohistology in formalin-fixed paraffin-embedded tissue sections. *Biochem Soc Trans* 1987;15:1179-80.
- Haines AMR, Larkin SE, Richardson AP, Stirling RW, Heyderman E. A novel hybridoma antibody (PASE/4LJ) to human prostatic acid phosphatase suitable for immunohistochemistry. *Br J Cancer* 1989;60:887-92.
- Nadji M, Tabei SZ, Castro A, Chu TM, Morales AR. Prostatic origin of tumors. An immunohistochemical study. *Am J Clin Path* 1980;73:735-9.
- Abrahamsson PA, Lilja H. Three predominant prostatic proteins. *Andrologia* 1990;22 Suppl 1:122-31.
- Allen SM. An enzyme linked immunosorbent assay (ELISA) for detection of seminal fluid using a monoclonal antibody to prostatic acid phosphatase. *J Immunoassay* 1995;16:297-308.
- Kämmerer U, Kapp M, Gassel IAM, Richter T, Tank C, Dietl J, et al. A new rapid immunohistochemical staining technique using the EnVision antibody complex. *J Histochem Cytochem* 2001;49:623-30.

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Agilent Technologies Singapore (International) Pte Ltd.