**Monoclonal Mouse
Anti-Human
Synaptophysin
Clone DAK-SYNAP**

Code M7315

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Synaptophysin, Clone DAK-SYNAP, is intended for use in immunohistochemistry (IHC). This antibody labels synaptophysin-expressing cells in normal and neoplastic formalin-fixed, paraffin-embedded (FFPE) tissue (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Synonyms for antigen	<p>Major synaptic vesicle protein p38, protein p38.</p>
Summary and explanation	<p>Synaptophysin is an acidic, homo-oligomeric, integral membrane glycoprotein (monomer = 38 kDa) originally isolated from rat brain (3) and belongs to a family of proteins with four transmembrane domains including synaptogyrin and synaptoporin (4). Despite being the most abundant synaptic vesicle membrane protein, the function remains enigmatic. Studies have shown numerous of diverse roles for synaptophysin such as to promote the formation of highly curved membranes in synaptic vessels, exocytosis, synapse formation, biogenesis and endocytosis of synaptic vessels as well as a role in learning and memory through regulation of synaptic transmission in neuronal circuits (4).</p> <p>Synaptophysin is localized to presynaptic vesicles of all neurons and in corresponding vesicles in neuroendocrine cells (5). Besides localization in neuroendocrine cells, synaptophysin is also found in non-neuroendocrine cells in adrenal cortical epithelia (3).</p> <p>Antibodies to synaptophysin have been shown to aid in the classification of neuroendocrine neoplasms, including neoplasms of epithelial type (6-8), nervous system neoplasms with neuronal differentiation (5, 6, 9) and adrenocortical neoplasms (10, 11).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant (containing fetal calf serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 0.015 mol/L sodium azide.</p> <p><u>Clone:</u> DAK-SYNAP. <u>Isotype:</u> IgG1, kappa.</p> <p><u>Mouse IgG concentration:</u> See label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
Immunogen	<p>A recombinant protein fragment corresponding to the C-terminal cytoplasmic domain of human synaptophysin.</p>
Specificity	<p>In Western blot analysis of human neuroblastoma cell lysate, the antibody labels a major band at 40 kDa corresponding to the expected molecular weight of synaptophysin.</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
Specimen preparation	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Tissue specimens should be cut into sections of approximately 4 µm.</p> <p><u>Pre-treatment:</u> Pre-treatment of formalin-fixed, paraffin-embedded tissue sections with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues with HIER using diluted EnVision FLEX Target Retrieval Solution, High pH (50x) (Code K8004). Deparaffinization, rehydration and epitope retrieval can be performed in Dako PT Link. For details, please refer to PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove slide rack from PT tank and immediately dip slides in jar/tank (e.g., PT Link Rinse Station (Code PT109)) containing diluted room temperature EnVision FLEX Wash Buffer (20x) (Code K8007). Leave slides in Wash Buffer for 1-5 minutes.</p> <p>The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended. After staining, the sections must be dehydrated, cleared and mounted using a permanent mounting method.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> The recommended dilution of Monoclonal Mouse Anti-Human Synaptophysin, Clone DAK-SYNAP, Code M7315, is 1:50. Dilute the antibody in Dako Antibody Diluent (Code S0809). Incubate pretreated tissue sections for 20 minutes at room temperature..</p> <p><u>Negative control:</u> The recommended negative control reagent is Dako Negative Control, Mouse IgG1 (Code X0931), diluted to the same IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, dilute these reagents immediately prior to use. Positive and negative controls should be run simultaneously with patient specimens.</p> <p><u>Visualization:</u> The recommended visualization system is EnVision FLEX, High pH (Code K8000/K8010) using a 20 minute incubation at room temperature. Follow the procedure enclosed with the selected visualization system(s).</p>

Automation: The antibody is well-suited for immunohistochemical staining using automated platforms, such as Dako Autostainer, Autostainer Plus and Autostainer Link as well as PT Link for pre-treatment.

Counterstaining: The recommended counterstain is EnVision FLEX Hematoxylin (Code K8008/K8018).

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include colon and the cells/structures should display reaction patterns as described for this tissue in the "Performance characteristics" section.

Staining interpretation

The cellular staining pattern is cytoplasmic.

Performance characteristics

Normal tissues: In colon, axons of the nerves in the tunica muscularis and in the lamina propria show a moderate to strong staining reaction, whereas the neuroendocrine cells in the epithelial surface show a weak to moderate staining reaction. Occasionally, cytoplasmic labeling of goblet cells in colon and small intestine can be observed.

Summary of Normal Tissue Reactivity (1).

Tissue Type (# tested)	Labeled Tissue Elements	Tissue Type (# tested)	Labeled Tissue Elements
Adrenal (3)	3/3 Adrenal cells (10-25%), cytoplasmic	Ovary (3)	1/3 Ganglion cells and stromal cells (10%) cytoplasmic
Bone marrow (3)	0/3	Pancreas (3)	3/3 Islet cells and ganglion cells (100%), cytoplasmic
Breast (3)	0/3	Pituitary (3)	2/3 Pituitary cells (100%) cytoplasmic
Cerebellum (2)	2/2 Ganglion and neuronal cells (100%), cytoplasmic	Placenta (3)	0/3
Cerebrum (3)	3/3 Ganglion and neuronal cells (100%), cytoplasmic	Prostate (3)	1/3 Stromal cells (50%), cytoplasmic
Cervix (3)	2/3 Ganglion cells (50%), cytoplasmic	Salivary gland (3)	3/3 Neuronal cells (100%), cytoplasmic
Colon (3)	3/3 Ganglion cells and neuroendocrine cells (100%), cytoplasmic	Skin (3)	2/3 Basal membrane cells (50%) diffuse
Endometrium (3)	0/3	Small intestine (3)	3/3 Neuroendocrine cells (100%), cytoplasmic
Esophagus (3)	3/3 Ganglion cells (100%), cytoplasmic 3/3 Fibroblasts (50%), cytoplasmic	Spinal cord (3)	3/3 Ganglion cells (100%), cytoplasmic
Fallopian tube (3)	0/3	Spleen (3)	0/3
Kidney (3)	0/3	Stomach (3)	3/3 Neuroendocrine cells (100%), cytoplasmic
Liver (3)	0/3	Testis (3)	0/3
Lung (3)	0/3	Thyroid (3)	0/3
Lymph node (3)	0/3	Tonsil (3)	0/3
Muscle, cardiac (3)	3/3 Fibroblasts (20%), cytoplasmic	Uterus (3)	0/3
Muscle, skeletal (3)	0/3	Ureter (3)	3/3 Stromal cells (5-10%) cytoplasmic
Nerve, peripheral (3)	3/3 Ganglion cells (10-100%), cytoplasmic	Urinary bladder (3)	0/3

Abnormal tissues: The antibody labeled 4/4 small cell lung carcinoma, 2/6 large cell lung carcinoma, 2/2 atypical lung carcinoid, 3/3 colon carcinoid, 5/5 small intestine carcinoid, 1/1 pancreas neuroendocrine neoplasia, 1/1 ganglioneuroma, 3/3 thyroid medullary carcinoma, 2/3 lung adenocarcinoma, 1/2 undifferentiated lung adenocarcinoma, 2/4 colon adenocarcinoma, 1/4 small intestine adenocarcinoma, 1/2 pancreas adenocarcinoma, and 0/3 thyroid, papillary carcinoma (2).

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References/ Bibliographie/ Literaturnachweise

- Dako in-house documentation D13338
- Dako in-house documentation D14505
- Wiedenmann B, Franke WW. Identification and localization of synaptophysin, an integral membrane glycoprotein of Mr 38,000 characteristic of presynaptic vesicles. Cell 1985;40:1017-28.
- Kwon SE, Chapman ER. Synaptophysin regulates the kinetics of synaptic vesicle endocytosis in central neurons. Neuron 2011; 70:847-85.
- Wiedenmann, B, Franke WW, Kuhn C, Moll R, Gould VE. Synaptophysin: a marker protein for neuroendocrine cells and neoplasms. Proc Natl Acad Sci U S A, 1986;83:3500-4.
- Gould VE, Wiedenmann B, Inchal L, Schwachheimer K, Dockhorn-Dvorniczak B, Radosevich B, et al. Synaptophysin expression in neuroendocrine neoplasms as determined by immunocytochemistry. Am J Pathol 1987;126:243-57
- Takei H, Asamura H, Maeshima A, Suzuki K, Kondo H, Niki T, et al. Large cell neuroendocrine carcinoma of the lung: a clinicopathologic study of eighty-seven cases. J Thorac Cardiovasc Surg 2002;124:285-92.
- Sorhaug S, Steinshamn S, Haaverstad R, Nordrum IS, Martinsen TC, Waldum HL. Expression of neuroendocrine markers in non-small cell lung cancer. APMIS 2007;115:152-63.
- Park SJ, Park CJ, Kim S, Jang S, Chi HS, Kim MJ, et al. Detection of bone marrow metastases of neuroblastoma with immunohistochemical staining of CD56, chromogranin A, and synaptophysin. Appl Immunohistochem Mol Morphol 2010;18:348-52.
- Li Q, Johansson H, Kjellmann M, Grimelius L. Neuroendocrine differentiation and nerves in human adrenal cortex and cortical lesions. APMIS 1998;106:807-17.
- Komminoth P, Roth J, Schroder S, Saremaslani P, Heitz P. Overlapping expression of immunohistochemical markers and synaptophysin mRNA in pheochromocytomas and adrenocortical carcinomas. Implications for the differential diagnosis of adrenal gland tumors. Lab Invest 1995;72:424-31.

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Agilent Technologies Singapore (International) Pte Ltd.