**Polyclonal Rabbit  
Anti-Human  
Thyroglobulin**  
Code A0251

**ENGLISH**

<b>Intended use</b>	<p>For in vitro diagnostic use.</p> <p>Polyclonal Rabbit Anti-Human Thyroglobulin, Code A0251, is intended for use in immunohistochemistry (IHC). The antibody labels thyroglobulin in thyroid tissue and is a useful aid for the classification of well differentiated thyroid carcinomas (1, 2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
<b>Summary and explanation</b>	<p>Thyroglobulin (Tg) is a glycoprotein with a predominant form as a 660 kDa homodimer. Tg, the precursor of thyroid hormones, is synthesized by thyrocytes and transported to the apical surface where it is secreted into the lumen of thyroid follicles and stored as the major component of colloid (&gt; 95%). A minor proportion of Tg is found as 330 kDa monomers or as tetramers. Reduction or degradation of 660 kDa or 330 kDa Tg molecules can lead to the formation of smaller polypeptides, some of which are present in trace amounts in the colloid. At the cell-colloid interface, post-translational modifications of Tg occur, which are characterized by coupling of tyrosyl residues with iodide, leading to the formation of thyroid hormone residues within the Tg molecule. Hormone release generally requires uptake of Tg from the colloid by thyrocytes and proteolytic cleavage along the lysosomal pathway (3).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure, Materials Required, Not Supplied, Storage, Specimen Preparation, Staining Procedure, Quality Control, Troubleshooting, Interpretation of Staining, General Limitations.</p>
<b>Reagent provided</b>	<p>Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L NaN<sub>3</sub>.</p> <p><u>Protein concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
<b>Immunogen</b>	Human thyroglobulin isolated from human thyroid glands.
<b>Specificity</b>	<p>The antibody labels human Tg, traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins.</p> <p>In crossed immunoelectrophoresis using 12.5 µL antibody per cm<sup>2</sup> gel area against 2 µL of human thyroid gland extract, only one precipitate corresponding to Tg appears. With 2 µL of human plasma, no precipitate appears. Staining: Coomassie Brilliant Blue.</p> <p>In human thyroid cell cultures, the antibody labels intracellular Tg (4).</p>
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For in vitro diagnostic use.</li> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> <li>Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
<b>Specimen preparation</b>	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Heat-induced epitope retrieval or proteolytic pre-treatment of deparaffinized tissues is not required. The tissue sections should not dry out during the immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> Labeling of thyroid cell preparations was obtained after fixation in 2% carbodiimide, Lison's "Gendre fluid", and 2 or 4% formalin, whereas Bouin, Carnoy A and B, formalin-calcium, Lillie's AAF, and especially 2% glutaraldehyde fixatives all gave a reduction in the amount of detectable intracellular Tg. No labeling was obtained after air-drying or treatment with methanol-H<sub>2</sub>O<sub>2</sub> or acetone (4). The user must validate the staining procedure.</p>
<b>Staining procedure</b>	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Polyclonal Rabbit Anti-Human Thyroglobulin, Code A0251, may be used at a dilution range of 1:4000-1:8000 when applied on formalin-fixed, paraffin-embedded sections of human thyroid tissue and using 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), Code X0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Visualization:</u> Dako EnVision+HRP kits, e.g. Code K4009, are recommended. Follow the procedure enclosed with the selected visualization kit.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p>
<b>Staining interpretation</b>	Cells labeled by the antibody display staining confined to the lumen of thyroid follicles and the apical surface of thyrocytes. In carcinomas the antibody also may display staining of the thyrocyte cytoplasm.
<b>Performance characteristics</b>	<p><u>Normal tissues:</u> The antibody labels Tg in thyroid follicles and at the apical surface of thyrocytes. Thyroglobulin is found in thyroid follicular cells and to a lesser extent, in interstitial tissue and thyroid stroma due to normal leakage.</p> <p><u>Abnormal tissues:</u> In 42 cases of primary thyroid carcinoma, labeling was observed with the antibody in 7/7 cases of papillary carcinoma (some cases were only stained along the apical surface, whereas others exhibited diffuse cytoplasmic staining) and in 16/16 cases of follicular carcinoma, regardless of the histological pattern, 10/14 cases of spindle and giant cell carcinomas were only weakly stained, and 0/5 cases of medullary carcinoma were labeled. All of 38 cases of carcinoma metastatic to the thyroid were not labeled, including those from kidney, lung,</p>

parathyroid, ovary, and breast, as well as the primary thyroid angiosarcomas and the soft tissue sarcomas of the neck that extended to the thyroid gland (1). In another study, the antibody was shown to label 10/10 cases of well differentiated thyroid carcinoma, five papillary and five follicular, whereas it labeled 0/53 cases of small cell anaplastic thyroid tumor, including 21 regarded as lymphomas, 13 as anaplastic carcinomas, and 19 as anaplastic tumors of uncertain histogenesis (2). The antibody labeled 0/16 cases of undifferentiated carcinoma, including 12 spindle-cell type and 4 mixed giant-cell and spindle-cell type, and 0/18 cases of malignant hemangiopericytoma (MHE) of the thyroid (6). In thyroid carcinoma of intermediate type, the antibody labeled more than 50% tumor cells in 9/18 cases, 25-50% tumor cells in 6/18 cases, less than 25% or scattered tumor cells in 3/18 cases (5).

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
**References / Bibliographie / Literaturnachweise**

1. Albores-Saavedra J, Nadji M, Civantos F, Morales AR. Thyroglobulin in carcinoma of the thyroid: An immunohistochemical study. Hum Pathol 1983;14:62-6.
2. Burt AD, Kerr DJ, Brown IL, Boyle P. Lymphoid and epithelial markers in small cell anaplastic thyroid tumors. J Clin Pathol 1985;38:893-6.
3. Marinò M, McLuskey RT. Role of thyroglobulin endocytic pathways in the control of thyroid hormone release. Am J Physiol Cell Physiol 2000;279:1295-1306.
4. Kayser L, Olsen BE, Høyer PE. Quantitative cytochemical demonstration of intracellular thyroglobulin in cultured human thyrocytes. Effects of fixatives, TSH and interleukin-1β. Histochem J 1991;23:235-40.
5. Ljungberg O, Nilsson PO. Intermediate thyroid carcinoma in humans and ultimobranchial tumors in bulls: A comparative morphological and immunohistochemical study. Endocr Pathol 1991;2:24-39.
6. Tötsch M, Dobler G, Feichtinger H, Sandbichler P, Ladumer D, Schmid KW. Malignant hemangioendothelioma of the thyroid. Am J Surg Pathol 1990;14:69-74.

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Agilent Technologies Singapore (International) Pte Ltd.